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
1. WHO. Report on the Burden of Endemic Health Care-Associated Infection Worldwide. [Accessed online on 14.03.2019]: [https://apps.who.int/iris/bitstream/handle/10665/80135/9789241501507\\_eng.pdf;jsessionid=86B33E75BBE594E6C2CBD998F1C88CA5?sequence=1](https://apps.who.int/iris/bitstream/handle/10665/80135/9789241501507_eng.pdf;jsessionid=86B33E75BBE594E6C2CBD998F1C88CA5?sequence=1)
2. WHO. Guidelines on Core Components of Infection Prevention and Control Programmes at the National and Acute Health Care Facility Level. [Accessed online on 14.03.2019]: [https://www.ncbi.nlm.nih.gov/books/NBK401773/pdf/Bookshelf\\_NBK401773.pdf](https://www.ncbi.nlm.nih.gov/books/NBK401773/pdf/Bookshelf_NBK401773.pdf)
3. WHO. Health Care-associated Infections: Fact Sheet. [Accessed online on 14.03.2019]: [https://www.who.int/gpsc/country\\_work/gpsc\\_ccisc\\_fact\\_sheet\\_en.pdf](https://www.who.int/gpsc/country_work/gpsc_ccisc_fact_sheet_en.pdf)
4. Al-Tawfiq JA and Tambyah PA. Healthcare associated infections (HAI) perspectives. J Infect Public Health. 2014;7(4):339-44.
5. Haque M, Sartelli M, McKimm J and Abu Bakar M. Health care associated infections - an overview. Infect Drug Resist. 2018;11: 2321-2333.
6. Agency for Healthcare Research and Quality. 2013 Annual Hospital-Acquired Condition Rate and Estimates of Cost Savings and Deaths Averted From 2010 to 2013. [Accessed online on 14.03.2019]: [https://www.ahrq.gov/sites/default/files/publications/files/hacrate\\_2013\\_0.pdf](https://www.ahrq.gov/sites/default/files/publications/files/hacrate_2013_0.pdf)
7. Dancer SJ, Stewart M, Coulombe C, et al. Surgical site infections linked to contaminated surgical instruments. J Hosp Infect. 2012;81(4):231-8.
8. Kirkland KB, Briggs JP, Trivette SL, et al. The impact of surgical-site infections in the 1990s: attributable mortality, excess length of hospitalization, and extra costs. Infect Control Hosp Epidemiol. 1999;20(11):725-30.
9. Teillant A, Gandra S, Barter D, et al. Potential burden of antibiotic resistance on surgery and cancer chemotherapy antibiotic prophylaxis in the USA: a literature review and modelling study. Lancet Infect Dis. 2015;15(12):1429-37.
10. Tosh PK, Disbot M, Duffy JM, et al. Outbreak of Pseudomonas aeruginosa surgical site infections after arthroscopic procedures: Texas, 2009. Infect Control Hosp Epidemiol. 2011;32(12):1179-86.
11. ECDC. Healthcare-associated Infections Acquired in Intensive Care Units. [Accessed online on 14.03.2019]: [https://ecdc.europa.eu/sites/portal/files/documents/AER-HCAI\\_ICU\\_3\\_0.pdf](https://ecdc.europa.eu/sites/portal/files/documents/AER-HCAI_ICU_3_0.pdf)
12. Timm D and Gonzales D. Effect of Sterilization on Microstructure and Function of Microsurgical Scissors. Surgical Services Management. 1997;3:47-49.
13. Rogers WJ. Steam and dry heat sterilization of biomaterials and medical devices. 2012:20-55.
14. Kovaleva J, Peters FT, van der Mei HC and Degener JE. Transmission of infection by flexible gastrointestinal endoscopy and bronchoscopy. Clin Microbiol Rev. 2013;26(2):231-54.
15. ASP Data on File. Why Go Green Now? The Benefits of Going EtO Free.
16. McGain F, Moore G and Black J. Steam sterilisation's energy and water footprint. Aust Health Rev. 2017;41(1):26-32.
17. McCreanor V and Graves N. An economic analysis of the benefits of sterilizing medical instruments in low-temperature systems instead of steam. Am J Infect Control. 2017;45(7):756-760.
18. ASP Data on File (2005). Methodist Hospital.
19. Shah SB and Bhargava AK. Recent advances in low temperature sterilization - Moving ahead from Cidex/ETO to OPA/Ozone: An update. Indian J Anaesth. 2017;61(10):855-857.
20. WHO. Concise International Chemical Assessment Document 54: ETHYLENE OXIDE. [Accessed online on 18.03.2019]: <https://www.who.int/ipcs/publications/cicad/en/cicad54.pdf>
21. WHO. Concise International Chemical Assessment Document 54: ETHYLENE OXIDE. [Accessed online on 18.03.2019]: <http://www.inchem.org/documents/cicads/cicads/cicad54.htm>
22. Ari Ş, Caca İ, Sahin A and Cingü AK. Toxic anterior segment syndrome subsequent to pediatric cataract surgery. Cutan Ocul Toxicol. 2012;31(1):53-57.
23. CDC. Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008. [Accessed online at on 05.04.2019]: [https://www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection\\_guidelines-H.pdf](https://www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection_guidelines-H.pdf)
24. U.S. Department of Energy FEMP. Water Efficiency Case Study - Veteran Affairs Hospital. 2007.
25. Rutala WA and Weber DJ. Disinfection and sterilization: An overview. Am J Infect Control. 2013;41(5):S2-S5.
26. McDonnell G and Burke P. Disinfection: is it time to reconsider Spaulding? J Hosp Infect. 2011;78(3):163-170.
27. Koo VS, O'Neill P and Elves A. Multidrug-resistant NDM-1 Klebsiella outbreak and infection control in endoscopic urology. BJU Int. 2012;110(11 Pt C):E922-6.
28. FDA. Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling. [Accessed online at on 18.03.2019]: <https://www.fda.gov/downloads/medicaldevices/deviceregulationandguidance/guidancedocuments/ucm253010.pdf>
29. FDA. Factors Affecting Quality of Reprocessing. [Accessed online on 18.03.2019]: <https://www.fda.gov/medicaldevices/productsandmedicalprocedures/reprocessingofreusablemedicaldevices/ucm454622.htm>
30. Dirlam Langlay AM, Ofstead CL, Mueller NJ, et al. Reported gastrointestinal endoscope reprocessing lapses: The tip of the iceberg. Am J Infect Control. 2013;41(12):1188-1194.
31. DiazGranados CA, Jones MY, Kongphet-Tran T, et al. Outbreak of Pseudomonas aeruginosa infection associated with contamination of a flexible bronchoscope. Infect Control Hosp Epidemiol. 2015;30(6):550-555.

asp.com

**ASP** Advanced Sterilization Products

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ADVANCED STERILIZATION PRODUCTS, INC.  
33 Technology Drive, Irvine CA 92618, USA

ASP, The Netherlands BV  
BIC 1, 5657 BX Eindhoven, The Netherlands

0123

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©ASP 2023. Tous droits réservés.
ASP France SAS